INTRODUCTION

Salvin[®] CytoSurg[®] Non-Resorbable PTFE membrane is available in several configurations and sizes. Provided that the Tyvek pouch is not compromised in any way, the Tyvek pouch will serve as an effective sterile barrier until the expiration date printed on the box. Salvin[®] CytoSurg[®] Non-Resorbable PTFE membranes should be stored in a cool, dry environment.

INDICATIONS

The Salvin[®] CytoSurg[®] Non-Resorbable PTFE Membrane is a temporarily implantable material (non-resorbable) for use as a space-making barrier in the treatment of periodontal defects.

DESCRIPTION

Salvin[®] CytoSurg[®] Non-Resorbable PTFE membrane is comprised of dense expanded polytetrafluoroethylene (ePTFE). ePTFE is recognized for being biologically inert and is tissue compatible. Through the use of a textured and microporous surface the ePTFE membrane is designed to provide a support scaffold for the regeneration of gingival soft tissue. Salvin[®] CytoSurg[®] Non-Resorbable PTFE membranes are designed to be a passive barrier and reduce the migration of epithelial and gingival tissue cells into the bony defect. The material is sufficiently firm to create and protect the defect space but supple enough to drape over and contour to the surrounding tissue.

MEMBRANE INSERTION

Carefully open the outer pouch and aseptically remove the inner pouch containing the membrane and place on to a sterile field. The sterile membrane can be removed from the second pouch during the surgical procedure. Handling of the membrane should only be done with sterile powder-free gloves or sterile atraumatic forceps. The membrane can be trimmed to size with sterile scissors allowing for adequate defect coverage. After trimming there should be no sharp corners or edges on the membrane. The membrane should be trimmed to extend 2-4mm beyond the defect margins and 1mm from adjacent, uninvolved teeth. If additional membrane stabilization after placement is needed, secure the membrane in place using a criss-cross or interrupted suture technique over top of the membrane, It is not recommended to suture through the membrane. The sutures are left in place for 10-14 days before being removed.

POST-OPERATIVE CARE

As with all oral surgical procedures a post-operative maintenance regiment is necessary for optimal healing. As prescribed by the clinician, the program may include rinsing with chlorhexidine in addition to gentle brushing with a soft bristle brush.

MEMBRANE REMOVAL

The Salvin[®] CytoSurg[®] Non-Resorbable PTFE membrane must be removed by 30 days maximum. If membrane is exposed, grasp with forceps and gently remove. If primary closure was achieved and membrane is not exposed a small incision will be necessary to remove membrane. Local anesthesia at the incision site will aid in patient comfort but is not necessary. With forceps, firmly grasp the membrane on the edge and gently work it back-and-forth if necessary, to aid in its removal. After membrane removal, re-epithelialization will occur within 14 to 21 days to complete the healing process. Depending on the patient, final bone maturation will take place over 6 to 12 months after placement.

CONTRAINDICATIONS

The Salvin[®] CytoSurg[®] Non-Resorbable PTFE membrane is a passive material and should not be used for load bearing applications. The Salvin[®] CytoSurg[®] Non-Resorbable PTFE membrane should not be considered for patients where general oral surgery is not advisable.


COMPLICATIONS

Possible complications associated with any oral surgery include infection, bone loss, abscess formation, flap slough, perforation, pain, and complications associated with anesthesia. Early removal of the membrane may be required depending upon clinician judgement of the severity and type of the complication.

 **CAUTIONS**

1. U.S.A. Federal Law restricts the sale, distribution or use of this device to, by or on the order of a licensed practitioner.
2. Do not use membrane if package has been damaged or opened prior to use.
3. Salvin[®] CytoSurg[®] Non-Resorbable PTFE membranes are single use only. Do not re-sterilize.
4. Salvin[®] CytoSurg[®] Non-Resorbable PTFE membranes are contained in sterile double Tyvek pouches and should be handled with sterile instruments in a sterile surgical field.
5. Salvin[®] CytoSurg[®] Non-Resorbable PTFE membranes should not be used in the presence of an active infection. Depending on the type of infection as judged by the clinician, early removal of the membrane may be required.
6. Salvin[®] CytoSurg[®] Non-Resorbable PTFE membranes are not designed to be placed in load bearing applications.

LABELING SYMBOLS

	CAUTION		USE BY:		CONSULT INSTRUCTIONS FOR USE	LOT	BATCH CODE
REF	CATALOG NUMBER		MANUFACTURER	STERILE EO STERILIZED USING ETHYLENE OXIDE			
	TEMPERATURE LIMITATION	R_x	PRESCRIPTION ONLY	 DO NOT RE-USE			
	DO NOT USE IF PACKAGE IS DAMAGED		DO NOT RESTERILIZE	UDI UNIQUE DEVICE IDENTIFIER			

SYMBOLS

 **Manufactured For:**
Salvin Dental Specialties, LLC
3450 Latrobe Drive • Charlotte, NC 28211
Toll Free: 800-535-6566 US & Canada
Website: www.salvin.com

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